



## **FREEHAND**

The practical garment holder

The FreeHand garment holder is a very practical aid for daily catheterisation. Pants or similar garments are pulled down and held securely by FreeHand, so that both hands are free to perform catheterisation.

FreeHand's tiny size enables it to fit into any bag. At home, FreeHand has it's own place using the included wallholder.



### **Directions for Use**

FreeHand, once fixed in place, makes sure that the undergarments are sufficiently isolated from the catheterisation.

Simply hook FreeHand into the opened pants and then hook it into the front edge of the chair or wheelchair. That way you have both hands free, and besides that, you avoid contact of the catheter with your clothes and thus reduce the risk of infections.

- 1 Application sketch for FreeHand
- 2 Application sketch for FreeHand when using the catheter system Libero with an attached urine pouch